

Honeywell Specialty Chemicals Seelze GmbH
Postfach 10 02 62, 30918 Seelze
Tel.: +49 (5137) 999-0, Fax: +49 (5137) 999-123
www.honeywell.com

Honeywell Specialty Chemicals Seelze GmbH • Postfach 10 02 62 • 30918 Seelze

Ihre Zeichen

Ihre Nachricht

Unsere Zeichen

Telefon: 05137 999-353

Seelze, 19.07.2019

Hydranal CoE

To Whomsoever It May Concern

This letter confirms that there is no change in the chemical composition of the Hydranal titrants when sold as Hydranal Smart Chemicals. The suffix “-SC”, representing “Smart Chemical”, added after the original catalog number of the titrants indicates that an RFID tag has been applied to the titrant bottle. This RFID tag enables certain data concerning the reagent to be read by the Mettler Toledo SmartChemical Reader but in no way alters the composition of the titrant itself.

Customers wishing to upgrade their Hydranal titrants to Hydranal Smart Chemicals can do so without any corresponding change in method.



Supervisor Hydranal Center of Excellence



Manager Quality Assurance